

# Metro Seal™

Tube Securement Device



The Hans Rudolph Metro Seal™ is a disposable, medical adhesive and thermoplastic elastomer, tube securement device. It holds an accessory tube in place on a patient's face, creating a seal between the mask, tube and patient's skin.

The Metro Seal™ can be used with any size V2™ oro-nasal mask when it is necessary to pass and secure an accessory tube (12 to 18 French) under the mask to the patient's nasal or oral orifices.

## Key Features

- Disposable, single patient use only.
- Skin safe adhesive and thermoplastic elastomer.
- Enables a seal between the mask, nasogastric tube and patient's skin.
- The accessory tubes involved may include but are not limited to nasogastric tubes and esophageal balloon catheters.



Find out more here

| Code   | Description | Box Qty |
|--------|-------------|---------|
| 669353 | Metro Seal™ | 25      |



eakinrespiratory.com

Distributed by  
Armstrong Medical Limited  
Wattstown Business Park,  
Newbridge Road, Coleraine,  
BT52 1BS, Northern Ireland.  
T +44 (0) 28 7035 6029  
coleraine@eakinhealthcare.com